



# 40ft



## STANDARD CLEAN ROOM CONTAINER



2.438 x 12.192 x 2.896 mm

|                                 |   |
|---------------------------------|---|
| EXTERNAL DIMENSIONS             | 2.438 x 12.192 x 2.896 mm   |
| INTERNAL DIMENSIONS (ISO8)      | Gowning Room: 2.130 x 1.200 x 2.200 mm<br>Laboratory Room: 2.130 x 10.060 x 2.200 mm  |
| INTERNAL DIMENSIONS (ISO6/ISO7) | Gowning Room: 2.130 x 920 x 2.200 mm<br>Laboratory Room: 2.130 x 9.780 x 2.200 mm   |
| CLEAN CATEGORY                  | ISO14644 - ISO8 – GMP GRADE D or ISO14644 - ISO6 / ISO7 – GMP GRADE C   |
| AIR CHANGES/HOUR                | ISO8 Gowning room: 10x Laboratory room: 15x<br>ISO6/ISO7 Gowning room: 15x Laboratory room: 60–90x  |
| KEY FEATURES                    | <ul style="list-style-type: none"> <li>Working area: 22 m2</li> <li>Suitable for 4–6 staff members</li> <li>Can be installed internally or externally</li> <li>Wind- and weatherproof intermodal container body</li> <li>Equipped with 2x cleanroom doors with EN179 emergency exit specs</li> <li>Interlocked entry system</li> <li>Complete clean HVAC system for air supply, filtration, and circulation</li> <li>Adjustable temperature and humidity control</li> <li>Manual positive pressure adjustment</li> <li>Touch-panel control board for HVAC system operation</li> <li>LED lighting throughout, 500 lux</li> <li>GMP grade vinyl floor</li> <li>Built-in wall sockets (AC and Ethernet)</li> <li>Comply with ISO 14644 and GMP standards</li> <li>Plug &amp; Play solution</li> <li>Electrical supply specified: 400V / 50Hz / 3~ / 3x16A</li> </ul> |
| MINIMUM HIRE PERIOD             | 16 weeks  |
| REQUEST A QUOTE                 | <a href="mailto:info@kleanlabs.com">info@kleanlabs.com</a>  |