

# 20ft


## CLEAN ROOM EXTENSION UNIT



EXTERNAL DIMENSIONS	2.438 x 6.058 x 2.896 mm
INTERNAL DIMENSIONS	2.130 x 5.280 x 2.200 mm
CLEAN CATEGORY	ISO6 – ISO7 or ISO8 (GMP C – GMP D) (can be stringer than main module)
AIR CHANGES/HOUR	Depending on cleanliness level (10x – 90x)
KEY FEATURES	<ul style="list-style-type: none"> <li>For connection with the main module (not standalone)</li> <li>Working area: 11 m<sup>2</sup></li> <li>Suitable for 2-3 staff members</li> <li>Can be installed internally or externally</li> <li>Wind- and weatherproof intermodal container body</li> <li>Equipped with 1x cleanroom door with EN179 emergency exit specs</li> <li>Door opening to the main module</li> <li>Complete clean HVAC system for air supply, filtration, and circulation</li> <li>Adjustable temperature and humidity control</li> <li>Manual positive pressure adjustment</li> <li>LED lighting throughout, 500 lux</li> <li>GMP grade vinyl floor</li> <li>Built-in wall sockets (AC and Ethernet)</li> <li>Comply with ISO 14644 and GMP standards</li> <li>Electrical supply specified: powered from main unit</li> </ul>
MINIMUM HIRE PERIOD	16 weeks
REQUEST A QUOTE	<a href="mailto:info@kleanlabs.com">info@kleanlabs.com</a>

# 40ft

## CLEAN ROOM EXTENSION UNIT



EXTERNAL DIMENSIONS	2.438 x 12.192 x 2.896 mm
INTERNAL DIMENSIONS	2.130 x 11.380 x 2.200 mm
CLEAN CATEGORY	ISO6 – ISO7 or ISO8 (GMP C – GMP D) (can be stringer than main module)
AIR CHANGES/HOUR	Depending on cleanliness level (10x – 90x)
KEY FEATURES	<ul style="list-style-type: none"> <li>For connection with the main module (not standalone)</li> <li>Working area: 25 m<sup>2</sup></li> <li>Suitable for 4–6 staff members</li> <li>Can be installed internally or externally</li> <li>Wind- and weatherproof intermodal container body</li> <li>Equipped with 1x cleanroom door with EN179 emergency exit specs</li> <li>Door opening to the main module</li> <li>Complete clean HVAC system for air supply, filtration, and circulation</li> <li>Adjustable temperature and humidity control</li> <li>Manual positive pressure adjustment</li> <li>LED lighting throughout, 500 lux</li> <li>GMP grade vinyl floor</li> <li>Built-in wall sockets (AC and Ethernet)</li> <li>Comply with ISO 14644 and GMP standards</li> <li>Electrical supply specified: powered from main unit</li> </ul>
MINIMUM HIRE PERIOD	16 weeks
REQUEST A QUOTE	<a href="mailto:info@kleanlabs.com">info@kleanlabs.com</a>